

# MCM 30"L EZ-Tranz Gold Patient Transfer Board | MCM 190



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 190 EZ-Tranz Gold Patient Transfer Board is a lightweight transfer solution weighing just 3 lbs, designed for safe lateral patient transfers between stretchers, beds, and surgical tables. Its conveyor-style design enables smooth, no-lift transfers while reducing caregiver strain.

The board features a soft, durable foam core for patient comfort and a gold antimicrobial cover with Ultra-Fresh protection for superior infection control. The surface allows for fast, easy cleaning with standard disinfectant wipes, supporting infection prevention in clinical environments.

Designed with integrated handles, the board is latex-free, MRI-safe, and supports patients up to 350 lbs. The 30" length accommodates standard transfer situations, making it ideal for operating rooms, imaging centers, and patient care floors.

## SPECIFICATIONS

- Length: 30"L
- Width: 15"L

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)