

MCM 48"L EZ-Tranz Gold Patient Transfer Board | MCM 191



Stock Photo - Options May Vary

FACTORY NEW

The MCM 191 EZ-Tranz Gold Patient Transfer Board is a 48-inch lateral transfer solution designed for safe, ergonomic patient movement across stretchers, beds, and surgical tables. It features a conveyor-style core that enables smooth, no-lift transfers while the soft, solid foam interior provides patient comfort.

The board includes a gold antimicrobial cover with Ultra-Fresh protection for enhanced infection control and cleans easily with standard hospital disinfectant wipes. It is MRI-safe, latex-free, and supports patients weighing up to 350 lbs, making it suitable for standard and bariatric transfer needs.

Integrated handles enable fast positioning and repositioning in clinical settings. The lightweight, durable construction reduces caregiver strain while the optional wall rack provides convenient storage and access in operating rooms, imaging suites, emergency departments, and other healthcare environments.

SPECIFICATIONS

- Length: 48"
- Width: 15"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)