

MCM Wall Rack for EZ-Tranz Gold Patient Transfer Board | MCM 199



Stock Photo - Options May Vary

FACTORY NEW

The MCM 199 Wall Rack provides efficient storage for EZ-Tranz Gold patient transfer boards in clinical settings. Constructed from stainless steel, it supports infection control efforts while withstanding daily use in demanding medical environments. The rack securely holds one board, keeping it off the floor and within arm's reach.

Designed to fit all EZ-Tranz Gold board sizes, this compact wall-mounted solution reduces clutter and improves workflow in operating rooms, imaging suites, and patient care areas. The rack does not include mounting hardware, allowing staff to select appropriate fasteners for specific wall surfaces such as drywall or concrete.

Once mounted, the MCM 199 offers a sturdy hold with a 25 lb weight capacity and is compatible across the EZ-Tranz line. Made in the USA and latex-free, the simple design ensures durability and ease of cleaning for hygienic storage and fast access to transfer boards.

SPECIFICATIONS

- Height: 2"
- Length: 2"
- Width: 6"
- Weight: Approx. 1 Pound

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)