

MCM Equalizing Device | MCM 350



Stock Photo - Options May Vary

FACTORY NEW

The MCM 350 Equalizing Device is designed for use when placing three or more digits in traction. It evenly distributes tension for consistent positioning and works with IV poles or extremity systems across diverse surgical environments.

Made with stainless steel hooks and a fully autoclavable chain, the MCM 350 delivers long-lasting performance and easy sterilization between uses. The device is ideal for orthopedic and hand procedures, ensuring reliable function while supporting proper technique and patient care.

Durable and easy to integrate into existing surgical setups, the MCM 350 is a valuable tool for any surgical setting requiring digit traction management.

FEATURES

- This device can be used with an I. V. pole or other extremity systems
- Stainless steel hooks and chain
- Fully Autoclavable

SPECIFICATIONS

- Stainless Steel
- Stainless Steel Hooks
- Stainless Steel Chain
- Fully Autoclavable

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)