

# MCM Foot Operated Mayo Stand | MCM 752



Stock Photo - Options May Vary

## FACTORY NEW

MCM 752: Foot Operated Mayo Stand. Mid Central Medical manufactures high-quality medical equipment to meet the demands of the American healthcare industry. The MCM 752 Mayo Stand offers a durable, user-friendly design ideal for hospitals, surgical centers, and outpatient facilities.

This stand includes solid brass inserts, which make height adjustment smooth and easy. Additionally, the unique "Safe Descent" feature lowers the tray gently, protecting patients and staff with controlled movement.

Furthermore, the removable 20-by-25-inch tray offers ample space for instruments and detaches easily for cleaning or replacement. The MCM 752 features a foot-pedal lock, enabling hands-free height adjustments to help maintain a sterile field during procedures.

## FEATURES

- Removable Tray
- Foot-Pedal Operated Lock
- Adjustable Height
- Low Profile Base Design
- Easily Slides Under Equipment

## SPECIFICATIONS

- Removable Tray Dimensions: 20" x 25"
- Adjustable Height Range: 36" to 62"
- Construction: All Welded #304 Stainless Steel
- Double Ball Bearing 2" Casters
- 2 Year Limited Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)