

MCM Hourglass Arm and Hand Surgery Table | MCM 325



Stock Photo - Options May Vary

FACTORY NEW

The MCM 325 is a carbon fiber arm and hand surgery table designed for operating room use. It features a single post adjustable leg and positive Posi-Lock clamps that secure firmly to the operating table, delivering unmatched stability. The extra-thick carbon fiber construction provides superior strength while maintaining a lightweight design at just 12 pounds.

The operating surface measures 15" x 32" and accommodates a wide variety of surgical procedures. The distinctive hourglass shape expands the table width to 24" at the ends, enhancing surgeon comfort and accessibility during procedures.

The table is fully radiolucent for clear imaging and includes a 2" thick pad. A built-in handle enables effortless carrying and serves as a convenient storage bracket. An optional Add-A-Rail can be attached to the end of the board to adapt multiple surgical accessories.

FEATURES

- Carbon Fiber Arm/Hand Table
- Lightweight
- Adjustable Single post leg
- Clamps positively to operating table with our unique Posi-lock clamps
- Made with extra-thick carbon fiber to provide superior strength when it is most essential
- Totally radiolucent for extremely clear images
- Hourglass shape expands the table width to 24" at the ends to enhance the surgeon's comfort
- The built-in handle allows effortless carrying and also acts as a convenient storage bracket, providing a place to hang the unit from when in storage

SPECIFICATIONS

- Weight: 12 lbs.
- The operating surface is 15" x 32"
- Width to 24" (at the ends)
- Leg adjustment range: 12"

**Brochure and images are a general overview of the product. Actual product options may vary.*