

MCM Steering Wheel I.V. Pole Accessory | MCM 230



Stock Photo - Options May Vary

FACTORY NEW

The MCM 230 is an IV pole steering wheel designed to improve mobility and control during patient care. It features rugged epoxy-coated steel construction that ensures long-term durability in fast-paced clinical environments. The wheel clamps securely to any 1-inch IV pole and tightens with two Allen screws for a strong, slip-resistant hold.

This steering wheel accessory helps healthcare professionals maneuver IV poles more easily and safely, especially in tight or crowded spaces. Its ergonomic design reduces strain and enhances precision when transporting IV equipment through hallways, patient rooms, or procedural areas.

Note: The IV pole is not included and must be purchased separately.

FEATURES

- Constructed with Rugged Epoxy Coated Steel
- Designed to Give Added Patient Stability when in Transport
- Clamps Anywhere on 1" Diameter Pole
- Tightens with Two Allen Screws to Give You a Positive Hold

SPECIFICATIONS

- Construction: Rugged Epoxy-coated Steel
- Application: Anywhere on 1" Diameter Pole
- Wheel Diameter: 16"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)