

MCM Tow Bar I.V. Pole Accessory | MCM 240



Stock Photo - Options May Vary

FACTORY NEW

MCM 240: IV Tow Bar. Mid Central Medical manufactures the MCM 240 IV Tow Bar to improve mobility and workflow in medical environments. Mid Central Medical proudly manufactures top-quality medical equipment to meet the demanding needs of healthcare professionals across the United States.

The MCM 240 features durable stainless steel construction for long-lasting performance in busy hospital and clinical settings. Moreover, the built-in clamp attaches to IV poles up to 1 1/4"—offering quick, tool-free setup and reliable stability.

Additionally, the MCM 240 easily adjusts from 20 to 30 inches in length, allowing flexible positioning based on patient transport needs.

FEATURES

- Stainless Steel Construction
- Integral Clamp
- Easily Clamps to IV Pole up to 1 1/4" Diameter
- Mounting Pin is 1/2" Diameter
- Adjustable Tow Bar

SPECIFICATIONS

- Stainless Steel Construction
- Clamps to IV Pole up to 1 1/4"
- Mounting Pin Diameter: 1/2"
- Adjusts from 20" to 30"
- Fits most Stretchers and Wheelchairs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)