

# MCM Mayo Tray Attachment | MCM 740



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 740 Mayo Tray Attachment by Mid Central Medical is a surgical accessory designed for operating room efficiency. It features stainless steel construction for durability and easy sterilization, with a removable tray design that allows for quick and thorough cleaning between procedures.

The attachment mounts securely using a 5/8" mounting post that fits into a standard Clark Socket (sold separately). Height adjustment up to 26 inches accommodates various surgical tasks and staff preferences.

The generous tray size provides ample workspace while maintaining unobstructed floor space, keeping the operating room organized and safe.

## FEATURES

- Unobstructed Floor Space in your O.R.
- Stainless steel construction for longevity
- The mounting post fits into a standard Clark socket (sold separately, see Option Tab)
- Generous tray size
- Height adjustable
- Removable tray for ease of cleaning

## SPECIFICATIONS

- Construction: Stainless steel
- Post Diameter: 5/8" mounting post
- Tray Dimensions: 12 5/8" x 19 1/8"
- Height Adjustable: Up to 26"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)