

MCM SS Height Adj. Instrument Table, 36 x 20 | MCM 530



Stock Photo - Options May Vary

FACTORY NEW

MCM 530 - SS Height Adjustable Instrument Table. Mid Central Medical is an American Company committed to Manufacturing the Highest Quality Medical Equipment to serve the Demands of America's Health Care Industry. The MCM530 Features: 20" Travel to accommodate many procedures.

The legs straddle the operating table for use during surgery. Furthermore, it also features a hand crank that adjusts the height, and its stainless steel fully welded construction has an extra reinforced top for added rigidity and durability.

The image may not be an exact representation of dimensions. Are you looking for a different size Instrument Table or more/less Shelves?

FEATURES

- 20" Travel to accommodate many procedures
- Legs straddle the operating table for use during surgery
- The hand crank adjusts the height
- Fully welded stainless steel construction
- Extra reinforced top for added rigidity and durability

SPECIFICATIONS

- Dimensions: 36"L x 20"W
- Leg Clearance: 28"
- Casters: 3" Locking Ball Bearing
- Crank: Manual Hand Crank
- Height Adjustment: 36" - 56"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)