

MCM SS Height Adj. Instrument Table, 36 x 24 | MCM 531



Stock Photo - Options May Vary

FACTORY NEW

MCM 531 - SS Height Adjustable Instrument Table. The MCM 531 offers strong, versatile support for surgeries. Made in the USA, it meets the high standards of healthcare workers. Additionally, the MCM 531 table adjusts up to 20 inches in height for easy use in many procedures.

Its legs fit around the operating table, saving space and improving access during surgery. The hand crank design makes adjusting the table's height simple and smooth, enabling quick changes without disrupting workflow. Medical staff appreciate this flexibility, which supports better ergonomics and patient care.

Furthermore, the MCM 531 is constructed from fully welded stainless steel. The MCM 531 features an extra-reinforced top for enhanced rigidity.

FEATURES

- 20" Travel to accommodate many procedures
- Legs straddle the operating table for use during surgery
- Hand crank adjusts the height
- Fully welded stainless steel construction
- Extra reinforced top for added rigidity and durability

SPECIFICATIONS

- Dimensions: 36"L x 24"W
- Leg Clearance: 28"
- Casters: 3" Locking Ball Bearing
- Crank: Manual Hand Crank
- Height Adjustment: 36" - 56"

**Brochure and images are a general overview of the product. Actual product options may vary.*