

MCM SS Height Adj. Instrument Table, 46 x 34 | MCM 533



Stock Photo - Options May Vary

FACTORY NEW

MCM 533 - Height Adjustable Instrument Table. The MCM 533 from Mid Central Medical delivers versatility and durability for modern surgical environments. As a proud American manufacturer, MCM produces high-quality medical equipment designed to meet the rigorous demands of healthcare professionals. The MCM 533 offers 20 inches of height adjustment to support various surgical procedures.

Its legs straddle the operating table for optimal positioning and access. A hand crank enables smooth, precise height changes, promoting workflow efficiency. Furthermore, made from fully welded stainless steel, it features an extra-reinforced top for added strength and durability.

Additionally, this sturdy design withstands heavy use in busy operating rooms, ensuring reliable performance.

FEATURES

- 20" Travel to accommodate many procedures
- Legs straddle the operating table for use during surgery
- Hand crank adjusts the height
- Fully welded stainless steel construction
- Extra reinforced top for added rigidity and durability

SPECIFICATIONS

- Dimensions: 46"L x 34"W
- Leg Clearance: 26"
- Casters: 3" Locking Ball Bearing
- Crank: Manual Hand Crank
- Height Adjustment: 36" - 56"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)