

# MCM Stainless Steel “Elite” Clark Socket | MCM 159



*Stock Photo - Options May Vary*

## FACTORY NEW

The MCM 159 is a stainless steel Clark socket designed for surgical table applications. It features a snap-on design that enables quick setup and removal, reducing preparation and cleanup time. The socket connects anywhere along the table rail, including between standoffs and other accessories, providing flexible positioning options.

This socket functions as two clamps in one, securely accepting both 3/4" diameter round posts and flat bar posts without requiring additional tools. Its all-#304 stainless steel construction provides exceptional durability and rust resistance, ensuring reliable performance in high-use surgical environments.

The MCM 159 delivers strength and versatility for demanding surgical settings. Its solid design and ease of use make it a reliable accessory that integrates well with modern surgical equipment and supports efficient workflow in operating rooms.

## FEATURES

- The new snap-on feature allows for quick setups and breakdowns
- Connects anywhere along the rail (even between standoffs and other accessories)
- The "Elite" is Two clamps in One
- Accepts flat blades and round posts

## SPECIFICATIONS

- Accepts flat blades and round posts up to 3/4" diameter
- Construction: Made with #304 Solid Stainless Steel

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)