

MCM SS Height Adj. Instrument Table, 48 x 24 | MCM 532



Stock Photo - Options May Vary

FACTORY NEW

MCM 532 - Height Adjustable Instrument Table. The MCM 532 by Mid Central Medical delivers versatile support for surgical teams across many specialties. This American-made table offers 20 inches of height travel, allowing precise positioning during diverse procedures. Its legs straddle the operating table, maximizing access and space efficiency in the operating room.

Medical professionals benefit from the enhanced flexibility and comfort this design provides. A hand crank controls the height adjustment, enabling smooth and accurate changes with minimal effort. This feature helps maintain surgical workflow without interruptions.

The crank's ergonomic design ensures easy operation for staff of all skill levels. As a result, the MCM 532 supports faster setup and improved efficiency.

FEATURES

- 20" Travel to accommodate many procedures
- Legs straddle the operating table for use during surgery
- The hand crank adjusts the height
- Fully welded stainless steel construction
- Extra reinforced top for added rigidity and durability

SPECIFICATIONS

- Dimensions: 48"l x 24"w
- Leg Clearance: 40"
- Casters: 3" Locking Ball Bearing
- Crank: Manual Hand Crank
- Height Adjustment: 36" - 56"

**Brochure and images are a general overview of the product. Actual product options may vary.*