

Stryker Stretcher Replacement Wheel Grey - New



Stock Photo - Options May Vary

FACTORY NEW

Stryker Stretcher Replacement Wheel Grey - New. Upgrade your stretcher with this brand-new 8-inch grey replacement wheel. Additionally, designed for durability, it ensures smooth movement across hospital floors.

The wheel comes with two dust caps for added protection and long-lasting performance. Compatible with multiple Stryker models, it replaces part numbers 0753-103-215 and 0715-002-025. Moreover, this non-OEM replacement delivers reliable performance without compromising quality or safety, making it ideal for medical facilities.

Furthermore, this replacement wheel fits: Stryker Prime 1105, 1115, and 1125 Stryker M Series models 1005, 1007, 1015, and 1025 Older Stryker models, including 0735, 0737, 0748, 1001, 1009, 1501, 1509, 1210, 1211, 1231, 1710, 1711, 1731, 721, 1061, 1068, 1069, 1079, and 1089. Crafted from durable materials, this replacement Stretcher Wheel withstands heavy loads and frequent use.

FEATURES

- 8" Stretcher Wheel (203mm)
- Grey
- Brand New
- Non-OEM Replacement
- Durable Materials

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)