

Stryker Stretcher Backrest Gas Spring | 0360-031-077 – New



Stock Photo – Options May Vary

FACTORY NEW

Stryker Stretcher Backrest Gas Spring: 420N. Improve backrest performance with this Stryker stretcher backrest gas spring, part number 0360-031-077. Designed for precise fowler adjustment, it delivers strong 420N support. Moreover, this brand new non-OEM replacement ensures consistent lifting power during daily use.

If your stretcher uses only one gas spring on the backrest, this strut likely fits your unit. As a direct replacement, it installs quickly and reduces downtime. Additionally, this gas spring fits all Stryker M Series stretchers.

Also fits the Stryker Prime series stretchers. Compatible models include Stryker 1005, 1007, 1015, and 1025. It also supports Stryker 1105, 1115, and 1125 models.

FEATURES

- Replaces Stryker P/N: 0360-031-077
- Used for fowler (backrest) adjustment
- Fits Stryker M Series and Prime Stretchers
- 420N
- Brand New

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)