

Stryker Rugged Ambulance Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Stryker Rugged Ambulance Stretcher is designed for emergency transportation with lightweight yet durable construction that simplifies operation and handling. One-handed operation of the backrest and side rails reduces user fatigue during patient transfers and transport.

The break-away head section reduces stretcher length by 62 inches, enabling navigation through tight spaces and confined areas. This design makes the Stryker Rugged one of the most maneuverable stretchers available for emergency services.

The stretcher features a 500 lb. weight capacity and comes fully refurbished with a new pad and warranty, ensuring reliable performance for ambulance and emergency transport applications.

FEATURES

- Transport at load height capability
- Positive action height adjustment
- Easy-to-use release handle design
- Color-coded controls
- Lightweight, rugged aluminum construction
- One-hand release, breakaway head section
- One-hand release, fold-down side rails
- One-hand release, infinite positioning, pneumatically assisted backrest

SPECIFICATIONS

- Overall Length: 62" - 80.5"
- Overall Width: 23"
- Weight: 81 lbs.
- Patient Surface Width: 26"
- Height Range: 13.5" - 37.5"
- Weight Capacity: 500 lbs.
- Trend/Rev Trend: +/- 18°
- Wheel Diameter: 6"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)