

MCM Stainless Steel Anesthesia Cart with Removable Tray | MCM 523-T



Stock Photo - Options May Vary

FACTORY NEW

The MCM 523-T: Anesthesia Cart w/ Removable Tray from Mid Central Medical is built for durability, organization, and efficiency in clinical settings. Constructed from high-quality stainless steel, it offers long-lasting performance and resistance to wear, making it ideal for demanding operating rooms and healthcare environments.

This cart has a large removable tray for easy access and moving important anesthesia tools. Additionally, it has four drawers to keep medicines, supplies, and small tools organized.

A lower shelf adds space for bigger equipment or bulky items, giving more room to store everything needed. The MCM 523-T supports a streamlined workflow by keeping critical supplies within reach during fast-paced procedures.

FEATURES

- Stainless steel construction for durability
- Large removable tray size
- 4 drawers
- Lower shelf for extra storage
- Casters, Two swivel casters

SPECIFICATIONS

- Construction: Stainless Steel
- Tray Dimensions: 16 1/4" x 21 1/4"
- Overall Dimensions: 16"deep x 20"wide x 34" High
- 4 drawers
- 1 Lower shelf

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)