

# MCM Stainless Steel Utility Table | MCM 521



Stock Photo - Options May Vary

## FACTORY NEW

MCM 521 - Stainless Steel Utility Table. The MCM 521 offers durability, function, and quality for modern medical settings. Mid Central Medical innovates its product line to meet changing customer needs.

Additionally, they use premium components to ensure lasting performance and reliability. The MCM 521 is made with strong 16-gauge stainless steel that resists rust, perfect for busy hospitals and clinics. It has two large drawers with smooth glides, so you can easily and safely access tools and supplies.

The smooth drawers help staff work faster during busy procedures. A three-sided guardrail surrounds the table's top, adding safety by keeping items securely in place during movement or use.

## FEATURES

- Stainless steel
- 3-sided guardrail
- 2 locking, 2 non-locking, 3" casters
- Two Drawers
- Drawers are mounted on stainless steel glides

## SPECIFICATIONS

- Construction: 16 gauge #304 Stainless Steel
- Dimensions: 16" x 20" x 34"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)