

# MCM Stainless Steel Utility Table | MCM 522



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 522 is a stainless steel utility table designed for hospitals, clinics, and healthcare facilities. It features durable 16-gauge #304 stainless steel construction for corrosion resistance and long-lasting strength. The table includes three large drawers mounted on smooth stainless steel glides for secure, organized storage of medical tools and supplies.

A three-sided guardrail surrounds the top surface to keep items secure and prevent spills or drops. The open-frame design allows for easy cleaning and maintenance of hygiene standards in clinical environments. Smooth-rolling casters enable quiet mobility for positioning the table wherever needed.

The MCM 522 combines robust construction with practical design to provide a dependable workspace with ample storage. It supports efficient workflows in busy medical environments while delivering durability and convenience.

## FEATURES

- Stainless steel
- 3-sided guardrail
- Three Drawers
- Drawers are mounted on stainless steel glides
- 2 locking, 2 non-locking, 3" casters

## SPECIFICATIONS

- Construction: 16 gauge #304 Stainless Steel
- Dimensions: 16" x 20" x 34"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*