

# MCM Stainless Steel Solution Stand | MCM 1000



Stock Photo - Options May Vary

## FACTORY NEW

MCM 1000 - Stainless Steel Solution Stand. The MCM 1000 from Mid Central Medical is built to meet the highest standards of durability and performance. As an American manufacturer, Mid Central Medical uses top-quality components to ensure reliability and consistency in every product.

This single basin stand features an H-brace design for added stability and uses durable 16-gauge #304 stainless steel. Additionally, the large 8½-quart capacity basin provides ample space for medical solutions and supplies. Equipped with 3-inch scuff-resistant casters, the stand moves quietly and smoothly across clinical floors, improving mobility and ease of use.

Its ergonomic height allows comfortable access, reducing strain during medical tasks.

## FEATURES

- Durable Construction
- Large Basin
- Rolls quietly along on 3" Scuff-Resistant Casters
- Rubber bumpers around the basin stops noise when rolling

## SPECIFICATIONS

- Construction: 16 Gauge #304 Stainless Steel
- Basin Capacity: 8-1/2 Qt
- 3" Scuff-Resistant Casters
- Dimensions: 15" x 18" x 34"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)