

MCM Stainless Steel Solution Stand | MCM 1002



Stock Photo - Options May Vary

FACTORY NEW

MCM 1002 - Stainless Steel Solution Stand. The MCM 1002 from Mid Central Medical delivers strength, reliability, and efficiency for daily clinical use. As a trusted American manufacturer, Mid Central Medical uses premium components to ensure consistent, high-quality performance that healthcare professionals can count on.

This double basin stand includes two large 8½-quart basins, offering ample room for solutions, instruments, and medical supplies. Additionally, the sturdy H-brace adds stability, and the 16-gauge stainless steel frame resists corrosion and wear in busy environments.

Additionally, the MCM 1002 includes 3-inch scuff-resistant casters, which let it roll quietly and smoothly across medical spaces.

FEATURES

- Durable Stainless Steel Construction
- Two Large basins
- Rolls quietly along on 3" Scuff-Resistant Casters
- H-Brace
- Rubber Bumpers around the basin stop noise when rolling

SPECIFICATIONS

- Dimensions: 15" x 33" x 34"
- Construction: 16 Gauge #304 Stainless Steel
- Basin Capacity: 8-1/2 Qt Each Basin

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)