

MCM Stainless Steel Solution Stand | MCM 1003



Stock Photo - Options May Vary

FACTORY NEW

MCM 1003 - Stainless Steel Solution Stand. The MCM 1003 Stainless Steel Solution Stand from Mid Central Medical offers durability, convenience, and mobility for everyday clinical use. As a trusted American manufacturer, Mid Central Medical builds each product using high-quality components to ensure long-lasting performance and reliability in demanding healthcare environments.

This double basin stand includes two large 8½-quart basins, providing ample space for holding medical solutions, instruments, or cleaning supplies. Additionally, the integrated lower shelf adds valuable storage space for extra tools, linens, or containers, helping staff stay organized in fast-paced settings.

Built with durable 16-gauge #304 stainless steel, the MCM 1003 resists rust and wear, even with frequent cleaning.

FEATURES

- Durable Stainless Steel Construction
- Two Large Basins
- Rolls quietly along on 3" Scuff-Resistant Casters
- With Shelf
- Rubber Bumpers around the basin stop noise when rolling

SPECIFICATIONS

- Dimensions: 15" x 33" x 34"
- Construction: 16 Gauge #304 Stainless Steel
- 3" Scuff-Resistant Casters
- Basin Capacity: 8-1/2 Qt (each basin)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)