

MCM Under Pad Mount Arm / Hand Surgery Table | MCM 340



Stock Photo - Options May Vary

FACTORY NEW

The MCM 340 is an under pad mount arm and hand surgery table designed for use during surgical procedures requiring C-arm access. It features a single post leg and a 2-inch thick pad, with a polycarbonate-covered opening measuring 14 inches wide by 24 inches long to accommodate imaging equipment.

The table is constructed from rugged stainless steel for durability and long-term use. Its design allows the table to slip under the operating table pad so patients do not feel its presence during procedures.

The table includes a convenient handle for easy carrying and storage. An optional Add-a-Rail can be attached to the end of the board to accommodate multiple accessories. The table ships complete with the pad.

FEATURES

- Very durable table
- Poly-Carbonate covered opening for C-arm access
- Rugged stainless steel construction for long-time use
- Table end designed to slip under operating table pad
- Thin enough that the patient will not feel it
- A convenient handle allows the table to be carried or stored effortlessly
- Optional Add-A-Rail can Easily be Added to the end of the board allowing you to adapt Multiple Accessories

SPECIFICATIONS

- Stainless Steel
- Opening Dimensions: 14"W x 24"L

**Brochure and images are a general overview of the product. Actual product options may vary.*