

MCM Under Pad Mount Arm and Hand Surgery Table with Tee Foot | MCM 341



Stock Photo - Options May Vary

FACTORY NEW

The MCM 341 is an under pad mount arm and hand surgery table designed for use during surgical procedures requiring C-arm imaging. It features a polycarbonate-covered opening measuring 14 inches wide by 24 inches long to provide clear access for precise imaging during procedures.

The table is constructed from rugged stainless steel for durability and long-term use. The design allows the table end to slip comfortably under the operating table pad so patients will not feel it during procedures. A convenient handle enables easy carrying and storage in busy clinical environments.

The MCM 341 includes a Tee Foot and a 2-inch thick pad. An optional Add-a-Rail can be attached to the board's end, allowing customization with multiple accessories to adapt for various surgical needs and improve workflow efficiency.

FEATURES

- Very durable table
- Polycarbonate-covered opening for C-arm access
- Rugged stainless steel construction for long-term use
- The table end is designed to slip under the operating table pad
- Thin enough that the patient will not feel it
- A convenient handle allows the table to be carried or stored effortlessly
- An optional Add-A-Rail (sold separately) can easily be added to the end of the board, allowing you to adapt multiple accessories

SPECIFICATIONS

- Stainless Steel
- Opening Dimensions: 14" x 24"

**Brochure and images are a general overview of the product. Actual product options may vary.*