

# MCM Wall Holder for 6" Wide Standard Armboards | MCM 445



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 445 is a wall-mounted holder designed for 6" wide standard armboards measuring 24 or 26 inches long. Constructed from durable, corrosion-resistant stainless steel, it securely holds armboards with pads attached while keeping them organized and accessible.

This wall holder saves valuable counter and floor space in busy medical environments by storing armboards out of the way when not in use. The smooth stainless steel construction protects armboards from scratches and damage, extending their lifespan while maintaining cleanliness and infection control standards.

Installation is quick and simple with included mounting hardware that fits most wall surfaces. The MCM 445 improves workflow and organization by keeping armboards within easy reach, making it ideal for hospitals, clinics, and other healthcare facilities.

## FEATURES

- Will hold armboards with the pads attached
- Keep your parts nearby and out of the way
- Keeps armboards handy and clean when not in use
- Will fit armboards 6" wide and 24 or 26 inches long

## SPECIFICATIONS

- Constructed of durable and easy to clean stainless steel
- Holds armboards 6" wide and 24 or 26 inches long

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)