

Medical Illumination Century Diagnostic Woods Hand Held Lamp | CWL-HHM



Century Diagnostic UV "Woods" Light
Stock Photo - Options May Vary

FACTORY NEW

Century Diagnostic UV "Woods" Light - Hand Held. The Century UV (Woods) Light offers exceptional value - in a full-featured light! Providing all the normal functions of a "Woods" light for specialties such as Dermatology or Ophthalmology, the Century also features two standard fluorescent bulbs for general illumination.

The traditional magnifier is augmented by an integral 6x "bifocal" lens for enhanced magnification. Additional Century models: MI also offers a variety of models to meet your needs.

Portable Floor Model with or without casters, Wall Mount, and Table/Shelf Clamp Mount - each featuring a quick-disconnect clamp for easy removal of light. Limited 3-Year Warranty. * Free shipping in the continental USA only.

FEATURES

- Exceptional value - in a full-featured light
- All the normal functions of a "Woods" light for specialties such as Dermatology or Ophthalmology
- Two standard fluorescent bulbs for general illumination
- The standard magnifier is augmented by an integral 6x "bifocal" lens for enhanced magnification
- Quick-disconnect clamp for easy removal of light

SPECIFICATIONS

- 6x "bifocal" lens
- 2x standard magnifier lens
- Two ultraviolet bulbs
- Two daylight fluorescent bulbs
- 6,000 hours rated bulb life

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)