

Medical Illumination MI 500 LED Exam/Diagnostic Single Ceiling Light | XLDE-SC



Stock Photo - Options May Vary

FACTORY NEW

Medical Illumination MI 500 LED Exam/Diagnostic Single Ceiling Light - XLDE-SC. The MI 500 light is part of the "cutting-edge" MI LED SERIES by Medical Illumination International. EFFICIENT: The MI 500's LED life averages 50,000 hours.

It uses significantly less electric consumption than standard halogen lights. COOL. The MI 500 LED Features: The light beam and head emit very low heat, making them part of our efficient LED solution.

Additionally, the MI 500 light boasts a brightness output of 50,000 lux (4,600 foot-candles) at a distance of 24", offering a whiter, more color-correct light. Despite its 6" diameter, its compact head design is powerful and offers excellent arm articulation.

FEATURES

- LED life averages 50,000 hours
- MI 500's light beam and head emit very low heat
- Light is output rated at 50,000 lux (4,600 foot candles) at a distance of 24.
- Whiter, a more color-correct light
- Its 6" diameter head design is small but powerful and has great arm articulation

SPECIFICATIONS

- Weight: Approximately 16.0 lbs
- Weight (Head): Approximately 3.5 lbs
- Weight (Arm): 12.5 lbs
- Arm: 57.0"
- Head Diameter: 6.25"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)