

Medical Illumination System Two Halogen Series | S2D-16



Stock Photo - Options May Vary

FACTORY NEW

The Medical Illumination System Two Halogen Series: S2D-16. (Two 16" Surgery Lights) incorporates proven halogen reliability and an unparalleled list of standard features, making it an industry standard.

Additionally, designed for outstanding performance, our System Two 16" halogen surgery light has features useful in various surgery and operating room settings. MI System Two Halogen Series S2D-16 Features: The Halogen Surgery Lights offer full focus control, 4-stage dimming, high color temperature, and various mounting options at a low cost.

Additionally, the redesigned System Two ceiling mounting and extension provide easier installation and improved stability, making it a more reliable and user-friendly solution for surgical lighting needs.

FEATURES

- up to 102,000 lux (9,500 foot candles) at 1 meter using a halogen bulb
- Over 4,000° Kelvin color temperature
- A fully adjustable light pattern
- Durable ceiling casting for optimum stability, balance, and function
- 360° rotational contacts at all critical points

SPECIFICATIONS

- Weight: Approximately 104.5 lbs
- Weight: Approximately 73.5 lbs
- Dimensions: Ceiling Casting: 17.0" diameter x 4.7" deep
- Rotations - Ceiling Mount/extension arm interface: Continuous
- Articulating arm vertical movement: +20 - 70 degrees

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)