

# Medical Illumination System Two Halogen Series | S2T-2016VC



Stock Photo - Options May Vary

## FACTORY NEW

Medical Illumination System Two Halogen Series: S2T-2016VC. Medical Illumination model number S2T-2016VC features System Two (1) 20" light, (1) 16" light, and (1) Video Camera.

With Medical Illumination's Trio, the widest range of capabilities exists. Furthermore, designed for outstanding performance, the S2T surgery light is useful in a variety of operating/surgery settings.

The MI System Two Halogen Series S2T-2016VC offers a three-surgery head model and combinations of two surgery lights/video cameras, lights/monitors, and more. The Trio provides complete surgery lighting with full video integration.

## FEATURES

- Many possibilities to meet the changing needs of the surgery room
- The flexibility of movement and central mounting at the same axis as the light allows for 300-degree rotation above the surgery table
- The arm is compatible with most medical-grade video monitors that have an S-Video input
- up to 102,000 lux (9,500 foot candles) at 1 meter using a halogen bulb
- Over 4,000<sup>o</sup> Kelvin color temperature

## SPECIFICATIONS

- Dimensions: Ceiling Casting: 17.0" diameter x 4.7" deep
- Rotations - Ceiling Mount/extension arm interface: Continuous
- Articulating arm vertical movement: +20 - 70 degrees
- Articulating arm/Yoke interface: Continuous
- Yoke/lamp head interface: +/- 150 degrees

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)