


Medical Illumination System Two Halogen Series | S2FM-20



Stock Photo - Options May Vary

FACTORY NEW

Medical Illumination System Two Halogen Series: S2FM-20. Medical Illumination's System Two Halogen Surgery Lights model number S2FM-20 (One 20" Surgery Light), incorporates proven halogen reliability as well as an unparalleled list of standard features in a Floor Model, making it an industry standard.

The System Two 20" halogen surgery light, S2FM-20, offers exceptional performance with full focus control, 4-stage dimming, high color temperature, and various mounting options at an affordable price. It provides excellent value for surgery and operating room settings.

Two light heads are available: the S2FM-20 with a 20" reflector (shown) and the S2FM-16 with a 16" reflector (not shown).

FEATURES

- Weight: Approximately 88.5 lbs
- Dimensions, Floor Base: 7.5" high x 28" long x 25" wide
- Dimensions, Upright: 2.5" diameter x 79" high
- Dimensions, Arm Assembly: 3.5" diameter (tapered) x 26.0" long
- Articulating arm vertical movement: +20 - 70 degrees

SPECIFICATIONS

- Weight, Floor Base with Upright: Approximately 88.5 lbs
- Arm Assembly: Approximately 16.3 lbs
- S2 Light Head: 17.5 lbs
- Dimensions, Floor Base: 7.5" high x 28" long x 25" wide
- Dimensions, Upright: 2.5" diameter X 79" high

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)