

# Medical Illumination System Two Halogen Series | S2S-20



Stock Photo - Options May Vary

## FACTORY NEW

The Medical Illumination System Two Halogen Series: S2S-20 incorporates proven halogen reliability and an unparalleled list of standard features, making it an industry standard. Additionally, the S2S-20 comes with a three (3) year manufacturer's limited warranty\*, ensuring long-lasting performance and peace of mind for users.

The System Two 20" halogen surgery light delivers excellent performance with full focus control, 4-stage dimming, and a 10,000-hour halogen bulb. Its multi-faceted reflector reduces shadows for flexible positioning.

Two light heads are available: the S2 with a 20" reflector and the S1 "satellite" with a 16" reflector. The redesigned ceiling mount improves installation and stability.

## FEATURES

- Proven halogen reliability
- Designed for outstanding performance
- Useful in a variety of surgery and operating room settings
- Full focus control
- 4 stage dimming

## SPECIFICATIONS

- Ceiling casting: 17.0" Diameter x 4.7" deep
- Ceiling rod: 2.5" Diameter x 9.0"-43.0" long (depending on ceiling height)
- Arm Assembly (to the center of the light head): 74"
- Light head assembly 23.0" Diameter x 9.5" deep
- Approximate Weight (Solo ceiling mount assembly): 73.5 lbs

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)