

# Medical Illumination System Two Halogen Series | S2D-20



Stock Photo - Options May Vary

## FACTORY NEW

The Medical Illumination System Two Halogen Series S2D-20 features two 20-inch halogen surgery lights with proven reliability and comprehensive standard features. The system delivers outstanding performance with full focus control, 4-stage dimming, high color temperature, and multiple mounting options. Two light head configurations are available: the S2 with a 20-inch reflector and the S1 satellite with a 16-inch reflector.

The S2D-20 provides 102,000 lux output at 1 meter distance with over 4,000° Kelvin color temperature and a fully adjustable light pattern. Electronic touch controls manage on/off and dimming functions, while a reserve bulb with auto-switchover ensures continuity in case of primary bulb failure. The system includes a removable sterilizable handle and a unique curved facial lens design.

The redesigned ceiling mount and extension arm provide easier installation and enhanced stability. The advanced suspension arm design delivers smooth, extensive coverage with 360° rotational contacts at all key points. Optional advanced video integration capability is available with video camera and monitoring options.

## FEATURES

- 102,000 lux (9,500 ft candles) output at 1 meter distance
- Over 4,000° Kelvin color temperature
- Reserve bulb with auto-switchover, in case of primary bulb failure
- Electronic touch controls for on/off and dimming
- Advanced suspension arm design, providing smooth extensive coverage

## SPECIFICATIONS

- Weight: Approximately 104.5 lbs
- Reflector: 20" diameter polished aluminum, faceted
- Central illuminance (adjustable): 2,000 -9,500 foot candles
- Focal length: 39"
- Operating temperature: 41° - 104°F

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)