

# Medical Illumination System Two Halogen Series | S2D-20MA



Stock Photo - Options May Vary

## FACTORY NEW

Medical Illumination System Two Halogen Series: S2D-20MA. The System Two M1 Monitor Arm is an ideal combination of the System Two series of surgery lights and advanced video imaging.

The result creates many possibilities to meet the changing needs of the surgery room. Note: Optional monitors are also available.

S2D-20MA Features Include its flexibility of movement and central mounting at the same axis as the light allows for 300-degree rotation above the surgery table. The arm is compatible with most medical-grade video monitors that have an S-Video input.

## FEATURES

- Many possibilities to meet the changing needs of the surgery room
- The flexibility of movement and central mounting at the same axis as the light allows for 300-degree rotation above the surgery table
- The arm is compatible with most medical-grade video monitors that have an S-Video input
- 102,000 lux (9,500 ft candles) output at 1 meter distance
- Over 4,000<sup>o</sup> Kelvin color temperature

## SPECIFICATIONS

- Please Note: Detailed specifications for the System Two Monitor Arm are currently under review.
- Weight: Approximately 73.5 lbs
- Dimensions: Ceiling Casting: 17.0" diameter x 4.7" deep
- Rotations - Ceiling Mount/extension arm interface: Continuous
- Articulating arm vertical movement: +20 - 70 degrees

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)