

# Medical Illumination System Two Halogen Series | S2D-20VC



*Stock Photo - Options May Vary*

## FACTORY NEW

Medical Illumination System Two Halogen Series: S2D-20VC. The System MI System is an ideal combination of the System Two series of surgery lights and advanced video imaging.

The result creates many possibilities to meet the changing needs of the surgery room. 8 input controllers are available, sold separately.

**S2D-20VC Features:** Its flexibility of movement and central mounting at the same axis as the light allows for 300-degree rotation above the surgery table. The Video Camera provides excellent video imaging capabilities, at an incredibly low price.

## FEATURES

- Many possibilities to meet the changing needs of the surgery room
- The flexibility of movement and central mounting at the same axis as the light allows for 300-degree rotation above the surgery table
- The arm is compatible with most medical-grade video monitors that have an S-Video input
- 102,000 lux (9,500 ft candles) output at 1 meter distance
- Over 4,000<sup>o</sup> Kelvin color temperature

## SPECIFICATIONS

- Weight: Approximately 104.5 lbs
- Dimensions, Ceiling Casting: 17.0" diameter x 4.7" deep
- Rotations - Ceiling Mount/extension arm interface: Continuous
- Articulating arm vertical movement: +20 - 70 degrees
- Articulating arm/Yoke interface: Continuous

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)