

Medical Illumination Centurion Excel Emergency Surgery Light w/Battery Backup | CH-FMB



Stock Photo - Options May Vary

FACTORY NEW

MI Centurion Excel Emergency Surgery Light w/Battery Backup. The Centurion Excel Emergency Surgery Light with Battery Backup (051513B), incorporates all of the fine features of the Centurion Excel Series, with an emergency backup system. The CH-FMB Emergency Surgery Light offers over 2 hours of uninterrupted light during power outages.

Its internal power control system ensures reliable lighting in adverse electrical conditions. The 16" multi-faceted reflector and color-correcting, heat-absorbing glass cylinder improve shadow reduction. The Centurion Excel delivers 64,500 lux (6000 foot-candles) at a 4,000 Kelvin color temperature for clear, precise illumination.

The CH-FMB is efficient: The 90-watt bulb is easily replaced after 2,000 hours of use.

FEATURES

- Provides over 2 hours of uninterrupted light output
- 16" diameter, multi-faceted reflector placed behind a color-correcting, heat-absorbing glass cylinder improves shadow reduction
- The Centurion Excel is output rated 64,500 lux (6000-foot candles) at 4,000 Kelvin color temperature
- Bulb life is 2,000 hour
- The 90-watt bulb is effortlessly replaced when needed

SPECIFICATIONS

- 6,000-foot candles (64,500 lux) output at 1 meter with 4,000° Kelvin color temperature
- Power supply circuitry, battery, and charging system are fully contained in a steel enclosure 19" W x 7.5" H
- The enclosure serves as the base for the lighting fixture and a 1.75" diameter upright tube for floor models
- Four Casters, two with locking mechanisms, allow easy mobility
- Operating voltage range: 110-125 VAC or 220-250 VAC

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)