

Stille Medstone3 CV Vascular Table

FACTORY NEW

The Stille Medstone3 CV Vascular Table is designed for high-performance imaging in endovascular and fluoroscopy procedures. It features smooth, free-floating movement for precise head-to-toe access, enhancing both speed and accuracy during positioning.



Stock Photo - Options May Vary

The table offers 36" longitudinal float and ± 4 " transverse movement, with a total imaging length of 64" and 32" for 3D imaging. Motorized height adjustment ranges from 26" to 42", providing versatility for various procedures while supporting patients up to 500 lbs.

Built for efficiency and reliability, the Medstone3 CV delivers high-quality imaging workflow performance with advanced positioning capabilities.

FEATURES

- Heavy duty steel base frame
- Four hospital grade locking casters
- Float control unit
- Mattress pad 2" thick with "Softcare™" technology
- One (1) hand control
- One (1) foot control
- One (1) catheter tray 24"
- Three (3) Patient safety straps

SPECIFICATIONS

- Tabletop dimensions 24-82"
- Longitudinal float 36"
- Transverse float ± 4 "
- Imaging length 64" head to toe
- 3D imaging length 32"
- Motorized elevation 26-42"
- Patient weight capacity 500lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)