

Biodex 840 Imaging Table – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Biodex 840 is designed for image-guided fluoroscopic procedures requiring stability, access, and precise positioning. It features a large radiolucent area free of cross members and a free-float X-Y tabletop that allows full fluoroscopic visualization and unobstructed C-Arm positioning. Motorized height control, isocentric lateral roll, and Trendelenburg motion provide comprehensive positioning capabilities.

The table supports a patient weight capacity of 500 lbs and includes a carbon fiber tabletop with stainless steel finish. Standard OR accessory rails are provided for secure patient positioning during procedures.

This completely refurbished unit comes with a new mattress and includes a 6-month warranty.

FEATURES

- Stable, Accessible, Precise style, Quiet, and Vibration-Free Positioning
- Carbon Fiber table top
- Stainless Steel Finish
- New Mattress
- Standard OR Accessory Rails
- Warranty Included

SPECIFICATIONS

- Tabletop Dimensions: 101" L x 24" W
- Overall Dimensions: 104" L x 26" W x 43.5" H (max height)
- Height: 33.5" to 43.5"
- Tabletop X Motion (head-to-toe float): 35"
- Tabletop Y Motion (side-to-side float): 10"
- Trendelenburg: 0° to 20°
- Isocentric Lateral Roll: 0° to ±15°
- Patient Capacity: 500 lb

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)