

Midmark 520 Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Midmark 520 is a refurbished medical stretcher designed for reliable patient transport and positioning in clinical settings. It features pneumatic-assist Fowler backrest adjustment up to 90 degrees, hydraulic foot-pedal height adjustment ranging from 26.5 to 36.5 inches, and Trendelenburg/reverse Trendelenburg positioning capabilities. The stretcher includes a central brake and steer system, fold-under side rails, and 8-inch casters for smooth maneuverability.

This stretcher accommodates patients up to 350 pounds and comes equipped with a new 4-inch foam pad. The unit has been completely refurbished to ensure reliable performance. Side rails measure 15 inches by 63.625 inches, providing secure patient containment during transport and care procedures.

FEATURES

- 350 Lbs Weight Capacity
- Hydraulic Foot Operated Height
- Fowler +/- (Gas Spring Assist)
- Trendelenburg/Reverse Trendelenburg
- Central Brake and Steer
- Fold Under Side Rails
- Completely Refurbished
- Includes New 4" Pad

SPECIFICATIONS

- Overall Dimensions: 80"L X 2
- Movement: Back Up/Down (Pneumatic Assist) 0° - 90°
- Trendelenburg/Reverse Trendelenburg: 12° / -7°
- Height Adjustable: 26½" - 36½"
- Weight Capacity: 350 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)