

Midmark 535 Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Cost-Effective, Completely Refurbished Midmark 535 Stretcher: Features Include: Pneumatic Assist Fowler (backrest) 90°, Hydraulic Height Adjustment (foot pedal on the foot end of the stretcher), Trend/Reverse Trend, Central Brake, and Fold-Down Side Rails. The stretcher has a 350 lbs Weight Capacity.

Midmark 535 Additional Features: The 535 has 8" Casters for easy maneuverability and control, 12" x 63 ¾" Side Rails, and each Midmark Stretcher comes with a New Pad and Warranty. Are you wondering about our Refurbishing Process?

Check out our blog, where we explain it all. 'Doing it Right: Explaining the Refurbishing Process.'

FEATURES

- Trendelenburg/Reverse Trendelenburg
- Back-Up/Down (Pneumatic Assist)
- Height Up/Down (Foot Pump at Foot End of Stretcher)
- Fold-Under Side Rails
- 4 Wheel Locking, Steering

SPECIFICATIONS

- Overall Dimensions: 84" x 33
- Back Up/Down (Pneumatic Assist); 0°-90°
- Height Up/Down: 25¾" - 35¾"
- Fold-Under Side Rails
- Brakes: 4 Wheel Locking

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)