

Midmark 540 Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Midmark 540 is a fully refurbished stretcher designed for use in operating rooms, emergency departments, post-operative care units, and patient transport. It features pneumatic-assist backrest adjustment up to 90 degrees, hydraulic foot-pedal height adjustment ranging from 24 to 33 inches, and Trendelenburg/reverse Trendelenburg positioning capabilities.

The stretcher accommodates patients up to 350 pounds and includes a central brake system with four-wheel locking and steering for secure positioning and smooth maneuverability. Eight-inch casters provide easy control in clinical environments, while fold-down side rails offer convenient access and safety.

Each unit comes with a new patient pad and includes a 6-month parts warranty. The stretcher is ready for immediate clinical use.

FEATURES

- Completely Refurbished
- Trendelenburg/Rev Trendelenburg
- Head Up/Down (Pneumatic Assist)
- Height Up/Down (Hydraulic)
- Fold Away Side Rails
- 4 Wheel Locking
- 350 Lbs Weight Capacity
- 6 Months Parts Warranty

SPECIFICATIONS

- Weight Capacity: 350 lbs
- Overall Length: 83"
- Overall Width: 29.5"
- Height Range - Low: 24"
- Height Range - High: 33"
- Fowler Backrest: 0° - 90°
- Patient Surface: 77" x 27"
- Casters: 8"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)