

Midmark 647 Power Podiatry Chair - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Midmark 647 Power Podiatry Chair - Refurbished. The refurbished Midmark 647 Power Podiatry Chair delivers comfort, precision, and efficiency, making it a smart choice for any podiatry practice. As a result, healthcare professionals continue to rely on Midmark's quality and performance to improve care and clinical efficiency.

Additionally, the Midmark 647 offers smooth electric height, back, and tilt adjustments, enabling providers to position patients easily and accurately. Its ergonomic design supports both patient comfort and provider access, helping make procedures more efficient and less physically demanding.

The Midmark 647 features: Low height (22.5"), which makes it easy and safe for patients to get on and off the table.

FEATURES

- Convenient, low height of 19" on standard base
- High height range of 40" on standard base
- Standard foot control makes positioning easy for the operator
- Foot control operates height, tilt, and back adjustments.
- Adjustable Foot extension

SPECIFICATIONS

- Max. Patient Weight Capacity 450 lb (204 kg)
- Minimum Height - 19" (48 cm) on Standard Base / 22.5" (57 cm) on Rotational Base
- Maximum Height - 40" (101 cm) on Standard Base / 43.5" (110 cm) on Rotational Base
- Width (Seat Section) - 29" (74 cm)
- Width (Foot Section) - 26" (66 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)