

Midmark Ritter 230 Power Procedure Chair / Table - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Midmark 230 is a universal power procedure chair designed for efficient patient care across multiple medical specialties. It features a height range from 22.5" to 40" with powered adjustments for backrest, leg rest, and tilt controlled via footswitch or hand control. The chair includes standard accessories such as an articulating headrest, side rails, paper holder, pull-out drain pan, and hideaway stirrups.

The foot/leg rest section locks in four proctology positions, accommodating diverse procedural needs. The chair is suitable for family practice, OB/GYN, internal medicine, dermatology, plastic surgery, colon and rectal procedures, and general surgery applications.

This refurbished unit meets electrical safety standards and supports a patient weight capacity of 450 lbs, making it a reliable choice for busy outpatient procedure environments.

FEATURES

- Table: Up/Down (22.5" - 40")
- Tilt (0° - 30°)
- Back: Up/Down (0° - 85°)
- Foot: Up/Down (0° - 90°)
- Foot Switch or Hand Control

SPECIFICATIONS

- Minimum Height: 22.5" (57.2 cm)
- Maximum Height: 40" (101.6 cm)
- Patient Weight Capacity: 450 lbs. (204.1 kg)
- Table Length (foot extended): 72" (185.4 cm)
- Width (seat): 28" (71.1 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)