

Midmark Ritter 630 Power Procedure Chair / Table – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Midmark 630 Power Procedure Chair – Refurbished. For decades, healthcare professionals have trusted Midmark products for their durability, quality, and reliable performance in delivering efficient patient care.

Furthermore, the Midmark 630 universal power procedure chair/ table is an excellent choice for multiple types of procedures. Midmark 630 Power Procedure Chair – Refurbished – Features: This chair boasts a shallow minimum 19" height.

This low entry and exit height makes it easy for the patient to safely and confidently get on and off the chair. Furthermore, this power chair comes with a footswitch or hand control for operating movements. In addition, the Midmark 630 comes standard with an accessory side rail, paper holder, and articulating headrest.

FEATURES

- Table: Up/Down
- Tilt
- Foot: Up/Down
- Back: Up/Down
- Adjustable Headrest

SPECIFICATIONS

- Minimum Height: 19" (48.26 cm)
- Maximum Height: 40" (101.5 cm)
- Table Length (foot extended): 72" (182.9 cm)
- Table Length (foot and headrest extended): 83" (210.8 cm)
- Width (seat): 28" (68.6 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*