

# Midmark Ritter 75L Power Procedure Chair / Table - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Midmark Ritter 75L is a power procedure chair/table designed for general exams, minor procedures, and OB/GYN applications. It features smooth electric movements with adjustable backrest, leg section, articulating headrest, and foot control for hands-free operation. The low minimum height of 22.5 inches facilitates safe patient transfer.

The chair includes an 8-way footswitch or hand control for operating all movements, accessory side rails, paper holder, and articulating headrest. Additional features include a pull-out drain pan, hideaway adjustable stirrups, and a lockable foot/leg rest section with four proctology chair positions.

This refurbished model is fully reupholstered and provides reliable performance for modern healthcare environments at a lower cost than new equipment.

## FEATURES

- Movement: Table Up/Down
- Tilt: Foot: Up/Down
- Back: Up/Down
- Adjustable Headrest
- Hideaway Adjustable Stirrups
- Pull-Out Cysto Pan
- Foot Switch or Hand Control
- Accessory Side Rails

## SPECIFICATIONS

- Minimum Height - 22.5" (57.15 cm)
- Maximum Height - 40" (101.6 cm)
- Table Length - 70.5" (179.1 cm) (extends to 81")
- Width (top) - 27" (68.6 cm)
- Tilt - Trendelenburg to 30 Degrees
- Foot or Hand Control - To power height, tilt, back foot/leg, Stop, auto-return
- Electrical - Domestic: 120 VAC, 60Hz, 12AMP, UL Listed, CUL Certified

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)