

Midmark Ritter 491 Power ENT Otolaryngology Chair – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Midmark Ritter 491 Power ENT Otolaryngology Chair – Refurbished. For decades, Health Care Professionals have come to depend on the quality, durability, and dependability of Midmark's products to enhance efficient patient care consistently. The Midmark 491 universal power procedure chair/ table is an excellent choice for multiple types of procedures.

Midmark Ritter 491 Power ENT Otolaryngology Chair Features: Low height ability of the chair, at only 22" Power controls are on both sides of the chair, making adjustments easy and convenient. Optional foot controls are also available. The back is operated at the touch of a button, and the height adjustment and the auto-return feature.

The base also rotates 330° for even more positioning flexibility.

FEATURES

- Low height ability
- Power controls are on both sides of the chair
- Optional foot control is also available
- Back, height adjustment, and auto-return feature is operated at the touch of a button
- Base rotates 330°

SPECIFICATIONS

- Minimum Height: 22" (55.9 cm)
- Maximum Height: 31" (78.7 cm)
- Width (Min): 24" (61 cm)
- Width (Max): 26.5" (67.3 cm)
- Rotation: 330°

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)