

Midmark 419 Power Exam Chair - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Midmark 419 Power Exam Chair - Refurbished. The refurbished Midmark 419 offers a dependable and cost-effective solution for today's medical professionals. Furthermore, for decades, healthcare providers have trusted Midmark products for their quality, durability, and long-lasting performance, helping to support efficient and effective patient care.

The Midmark 419 is a versatile power examination chair/table, ideal for a variety of procedures including general exams, dermatology, minor surgeries, and OB/GYN. Additionally, it features smooth electric adjustments, allowing for easy positioning to improve workflow and patient comfort.

Furthermore, the 419's compact design and adjustable positions make it ideal for practices needing to save space without losing functionality.

FEATURES

- Low height ability
- Footswitch controls make adjustments easy
- Back, height adjustment, and auto-return feature is operated at the touch of a button
- Movements and Features: Table Up and Down, Back-Up and Down, Tilt, Auto Return, and Arm Rests

SPECIFICATIONS

- Minimum Height: 22.4" (56.9 cm)
- Maximum Height: 40.4" (102.6 cm)
- Width: 23" (54.4 cm)
- Length: 65" (165.1 cm)
- with headrest extended adds 15" (38.1 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)