

Midmark 625 Power Exam Table - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Midmark 625 Power Exam Table - Refurbished. The Midmark 625 Barrier-Free Exam Table offers versatile features to support a wide range of clinical needs. The design creates a more efficient and effective caregiving and patient experience.

Midmark 625 Exam Table - Features: This exam table has a minimum height of 18", and allows for easy entry and exit by even the most challenged patients. The table also raises to a maximum of 37".

The Midmark 625 accommodates patients up to 650 lbs, preserving their dignity while you provide them with the best care. This exam table comes with a footswitch or hand control for operating movements.

FEATURES

- Table: Up/Down (18" - 37")
- Backrest: Up/Down
- Pelvic Tilt (7°)
- 650 Lbs Weight Capacity
- Pull Out Foot Section

SPECIFICATIONS

- Minimum Height: 22.5" (57.2 cm)
- Maximum Height: 37" (94 cm)
- Table Length (footrest extended): 75.5" (190 cm)
- Table Length (footrest stowed): 59.5" (151 cm)
- Patient Weight Capacity: 650 lbs. (295 kg)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)