

Midmark 223 Power Exam Table - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Midmark 223 is a barrier-free power exam table designed with safety and accessibility as primary considerations. It features powered height adjustment ranging from 18" to 37" and accommodates patients up to 400 lbs, making it suitable for a wide range of patient needs. The table includes powered backrest adjustment, pull-out leg section, drain pan, and hideaway stirrups.

The 223 offers multiple control options with foot or hand controls for height and backrest movements. Additional features include foot extension shelves to expand table length, a paper roll holder, and a side electrical receptacle. The table has been completely refurbished and reupholstered.

Key movements include powered height adjustment, powered backrest control, and a pull-out leg section. The design prioritizes ease of patient transfer and reduces injury risk compared to standard box-style exam tables.

FEATURES

- Height Adjustable (18" - 37")
- 400 lb. Weight Capacity
- Powered Backrest: Up/Down
- Hideaway stirrups that lock in one of four lateral positions
- Foot or Hand Control

SPECIFICATIONS

- Minimum Height: 18" (45.7 cm)
- Maximum Height: 37" (94 cm)
- Patient Weight Capacity: 400 lbs. (181.4 kg)
- Overall Width: 30.5" (77.5 cm)
- Back Adjustment: 0° to 80°

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)