

Midmark 222 Power Exam Table - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Midmark 222 Power Exam Table - Refurbished. The Midmark 222 Barrier-Free Exam Table prioritizes safety with its purposeful, patient-focused design.

Designed to reduce the likelihood of injury to patients. Additionally, the 222 will help those who might have difficulty moving on and off the industry-standard box-style exam table.

The Midmark 222 Exam Table is helping to create a new industry standard with powered height adjustment from just 18" to 37" and a patient weight capacity of up to 400 lbs. Furthermore, this power exam table comes with a footswitch for operating movements.

FEATURES

- Height Adjustable (18" - 37")
- Lift Assist Backrest: Up/Down
- Pelvic Tilt (7°) (optional feature)
- 400 lb. Weight Capacity
- Hideaway stirrups that lock in one of four lateral positions

SPECIFICATIONS

- Minimum Height: 18" (45.7 cm)
- Maximum Height: 37" (94 cm)
- Lift assist Back Adjustment: 0° to 80°
- Overall Width: 30.5" (77.5 cm)
- Table Width 28.5" (72.4cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)